



SIZING GUIDE HADDENHAM VENEX, STAR COTTON AND DOKTUS

This form should NOT be used for ordering Haddenham Custom Garments. For Custom Upper Limb Garments please use form F0032. Available at: www.hadhealth.com/forms

LENGTH



Take the length measurement on the inside of the arm. Measure from the wrist crease to the underarm (axilla).

c - g cm

Venex

LENGTHS (cm)

	Regular	Long
c - g	38 - 44	44 - 50


STAR Cotton

	Regular	Long
c - g	38 - 44	44 - 50

Doktus

	Regular	Long
c - g	40 - 44	44 - 50

CIRCUMFERENCES




a	Palm at base of fingers	<input type="text"/> cm
b	Palm at fold of thumb	<input type="text"/> cm
c	Wrist crease	<input type="text"/> cm
d	Widest point of forearm	<input type="text"/> cm
e	Elbow fold	<input type="text"/> cm
g	Under arm (axilla)	<input type="text"/> cm

CIRCUMFERENCES (cm)

	XSmall	Small	Medium	Large
a	17 - 18	18 - 20	20 - 22	22 - 24
b	18 - 19	19 - 21	21 - 23	24 - 26
c	17 - 18	18 - 20	20 - 22	22 - 24
d	23 - 24	24 - 26	26 - 28	29 - 31
e	23 - 25	25 - 27	28 - 30	30 - 32
g	28 - 32	32 - 36	36 - 40	40 - 44
g Extra Wide	33 - 36	36 - 40	40 - 44	44 - 48

IF PATIENT DOES NOT FIT INTO A SIZED GARMENT THEN A HADDENHAM CUSTOM GARMENT MAY BE MORE APPROPRIATE

SIZING GUIDE MICROFINE GLOVE



	XSmall	Small	Medium	Large	
a	Palm at base of fingers	17 - 18	18 - 20	20 - 22	22 - 24
b	Palm at fold of thumb	18 - 19	19 - 21	21 - 23	24 - 26